


## Looking after your Larval Therapy


### BioBag daily care

- Do not immerse in water. Do not occlude
- Avoid sustained, direct pressure as this may occlude the larvae. Short periods for the purposes of mobilisation are permissible
- Daily change of the secondary dressings is recommended and when strikethrough is present
- Re-apply barrier where necessary to the peri-wound area (**fig.1**)
- Check larvae are viable at secondary dressing changes – movement of larvae and presence of dark red exudate indicate the larvae are alive (**fig.2**)
- Ensure damp gauze and an absorbent pad are replaced on top of the BioBag at each secondary dressing change (**fig.3**)
- Ensure that all outer/ dressings are not occlusive and are permeable to the air (**fig.4**)
- After 72 hours, reassess wound to decide on further treatment. If a further NEW larval treatment is required, schedule a new order
- If debridement is near completion and no further NEW larval treatment is required, plan follow on care/dressings as per specialist instructions or local formulary
- On removal, double bag and treat as clinical waste in line with your local Grade A Clinical Waste Disposal Protocol.


### 4 Day Treatment Cycle




**Fig.1** Re-apply barrier cream




**Fig.2** Check larvae are viable, reposition BioBag if necessary



**Fig.3** Place a moistened gauze swab over the BioBag



**Fig.4** Secure the moistened swab and absorbent pad



## How to store and order

### BioBag sizes and codes

	Code	Description
BioBag	BB50	2.5x4cm
	BB100	5x4cm
	BB200	5x6cm
	BB300	12x6cm
	BB400	10x10cm

### Ordering Larvae

**Orders received by us before 2pm will qualify for inclusive next day delivery, or a future planned date of your choosing.**

Please allow time for your own internal procurement/pharmacy to process the order.

Telephone:

**0345 230 1810**

E-mail : [orders@biomonde.com](mailto:orders@biomonde.com)

Fax : 01656 668 047

### Office Hours

Monday to Friday 8:30am – 5:00pm

For assistance outside working hours please call our **Clinical Helpline: 0345 230 6806**

### Storage

- Keep in transit containers
- Do not freeze, do not store above 25°C
- Must be applied by expiry date; usually the day after delivery. For optimal results apply on day of delivery.

[www.biomonde.com](http://www.biomonde.com)