


Getting Started

1. Switch power on (rear of CardioQ-ODM).
2. Connect probe to Patient Interface Cable (PIC).
3. Press **New patient**.
4. Enter Patient ID number or press **Auto number**.
5. Select **Male** or **Female**.
6. Use the large control knob to dial in patient age. Press the control knob to enter. Repeat process for height and weight.
7. Refer to operating handbook if patient data is outside nomogram limits. Check and press **Accept data**.
8. Apply water-based lubricant to lower part of probe and insert into oesophagus.
9. For oral placement advance probe until incisors are at the second depth marker. When using nasal placement advance probe gently until nasal septum is at the third depth marker (nearest connector).

Never use excessive force to insert the probe as this may harm the patient.







Using the CardioQ-ODM - Getting Started

Locating the Descending Aortic Waveform


10. When locating the CardioQ-ODM signal adjust the volume knob as required.

11. Adjust probe depth to locate the descending aortic signal and then rotate to optimise the signal.




Using the CardioQ-ODM - Locating the Descending Aortic Waveform

Optimise the Waveform




12. Press **Peak velocity display**.

Adjust the probe to find highest blue line (peak) and the sharpest audible pitch to obtain the best signal quality.



13. Activate **Auto gain**.

Yellow line confirms auto gain activation.




14. Monitoring begins.

Green line and white arrows confirm start of monitoring.

Set vertical range to suit height of waveform. Press **Focus** and then press **Range** to select desired vertical range.

Using the CardioQ-ODM - Optimise the Waveform

The CardioQ-ODM Waveform



The **green line** indicates the velocity/time envelope which the monitor uses to make calculations. The **white arrows** indicate time and velocity values used for CardioQ-ODM calculations.

The Stroke Distance (SD) is the area under the waveform and is the basic measured parameter upon which calculations of Stroke Volume (SV) and all other Cardiac Output (CO) and indexed measurements are made. Stroke Volume is the parameter of choice for fluid management protocols, however Cardiac Index (CI) and Stroke Volume Index (SVI) can also be utilised.

The waveform base, (flow time) depends on heart rate, left ventricular filling and afterload. The flow time corrected to a heart rate of 60bpm (FTc) is inversely correlated with the systemic vascular resistance (SVR).

The most common cause of a short FTc (<330 ms) is hypovolaemia. If a short FTc (<330 ms) does not increase after an appropriate fluid challenge, other causes of vasoconstriction, (e.g. excess vasopressors, cold temperature, or obstructed circulation such as pulmonary embolus) should be considered. A long FTc (>360 ms) is seen in conditions associated with a low SVR e.g. sepsis, pregnancy, vasodilators and some anaesthetic agents.

Peak Velocity (PV) and Mean Acceleration (MA) are markers of left ventricular contractility. Changes in afterload will also affect the PV and MA - both decrease with an increase in afterload, and vice versa. Changes in preload predominantly affect the FTc and only affect PV and MA at extremes.

Using the CardioQ-ODM - The CardioQ-ODM Waveform

Basic Setup & Haemodynamic Parameters*

	SV ♦ Stroke Volume	Blood volume ejected during each systolic phase (ml).
	SD ♦ Stroke Distance	Distance a column of blood moves through the descending thoracic aorta during each systolic phase (cm).
	FTc Flow Time Corrected	Systolic flow time corrected for heart rate (ms).
	PV Peak Velocity	Peak velocity of blood flow in systolic phase (cm/s).
	SVI ♦ Stroke Volume Index	Stroke Volume normalised for body surface area (l/min/m²).
	SD ♦ Stroke Distance	Distance a column of blood moves through the descending thoracic aorta during each systolic phase (cm).
	FTc Flow Time Corrected	Systolic flow time corrected for heart rate (ms).
	PV Peak Velocity	Peak velocity of blood flow in systolic phase (cm/s).

Other parameters

CO Cardiac Output	Litres of blood pumped per minute (l/min).
CI Cardiac Index	Cardiac output normalised for body surface area (l/min/m²).
MD Minute Distance	Distance a column of blood moves through the descending thoracic aorta per minute (cm); MD = SD x HR; linear cardiac output.
HR Heart Rate	Beats per minute (bpm).
MA Mean Acceleration	Average acceleration of blood from start of systole to detected peak velocity (cm/s²).
SVR Systemic Vascular Resistance	The resistance that the left heart pumps against; measure of left ventricular afterload; note: external blood pressure data required to calculate SVR (dyn.s/cm⁵).
SVRI Systemic Vascular Resistance Index	Systemic vascular resistance normalised for body surface area (dyn.s/cm⁵/m²).

♦ Fluid management algorithm protocol parameters ♦ Alternative fluid management algorithm protocol parameters

* Refer to operating handbook for additional haemodynamic parameters provided by the CardioQ-ODM.

Using the CardioQ-ODM - Basic Setup & Haemodynamic Parameters




CardioQ-ODM Quick Reference Guide Intensive Care Application - Interpreting Results

This document is to be used by experienced clinical staff and is not intended to replace the operating handbook

Fluid Management

Typically, optimisation is achieved through the use of a fluid management algorithm. Stroke Volume (SV) or Stroke Distance (SD) responses to fluid challenges may help guide further interventions. Such algorithms have been utilised routinely in outcome studies with Deltex Medical Oesophageal Doppler Monitors.



Treatment algorithm suggested by Prof. M. Singer, University College Hospital London.

Interpreting Results - Fluid Management

Typical Parameter Values*

*These values should not be confused with a physiological target.

Flow Time Corrected (FTc)

330 - 360 milliseconds

Note – The effects of vasodilating drugs may elevate the FTc.

1. Singer, M. Oesophageal Doppler monitoring of aortic blood flow: beat by beat cardiac output monitoring. **International Anesthesiology Clinics** 1993; Vol.31(3):99-125.

2. Gardin, JM, Davidson, DM, Rohan, MK, et al. Relationship between age, body size, gender and blood pressure and Doppler flow measurements in the aorta and pulmonary artery. **Am Heart J** 1987; 113:101-109.


Peak Velocity (PV)^{1,2}

20 years	90 - 120 cm/s
30 years	85 - 115 cm/s
40 years	80 - 110 cm/s
50 years	70 - 100 cm/s
60 years	60 - 90 cm/s
70 years	50 - 80 cm/s
80 years	40 - 70 cm/s
90 years	30 - 60 cm/s

Extrapolated values are in plain text.


Interpreting Results - Typical Parameter Values

Intervention




Possible hypovolaemia.
Reduced SV/SD, decreased FTc.

FLUID




Positive fluid response.
SV/SD increases by more than 10%.




Left ventricular failure.
Reduced SV/SD and PV. Rounded waveform apex.

INOTROPE




After inotrope.
Increased SV/SD and PV.



Probable high SVR/afterload.
Reduced SV/SD, PV and FTc.

VASODILATE



Reduced SVR/afterload.
Increased SV/SD, PV and FTc.

Interpreting Results - Intervention

Additional Features

Number of Cycles per Calculation

To change cycle setting: While in the Run Mode, press **Cycles**. Rotate the large control knob to make cycle selection, and press the control knob to finish. Increasing cycle setting may aid parameter averaging on patients with an irregular rhythm or a respiratory swing, while decreasing cycle setting may be useful for monitoring during diathermy.

Storing a Waveform/Snap Function

While in a Run screen, press **Freeze**. Rotate the large control knob to place the desired section within the red Snap Window box. Press **Take snap**. The snap is then displayed in the split screen. To return to full screen press **Home** then **Full screen** and then **Run**. Record and view up to eight recorded waveform images.

Setting the Signal Filter

In Run Mode press **Focus**, then press **Filter** to activate the filter. Press **Filter** a second time to deactivate. Activate only to reduce artefact from excess heart valve or wall motion noise.

Refer to the CardioQ-ODM operating handbook for additional information.



Deltex™ and CardioQ-ODM™ are trademarks of Deltex Medical.
© Deltex Medical 2010. 9051-5514 - Issue 2
www.deltexmedical.com
CE 0086

Interpreting Results - Additional Features



CardioQ-ODM Quick Reference Guide Intensive Care Application - Getting Started

This document is to be used by experienced clinical staff and is not intended to replace the operating handbook