



## **Direct Current Cardioversion (DCCV)** (See overleaf for algorithm)

1. **Sedation** or **general analgesia** for conscious patients
2. With the defibrillator in manual mode, attach **3 lead** and defibrillation pads

3.  Press the **soft key “Sync On/Off”** and observe sense markers ↓

4.  If necessary, press  to select lead with most sense markers ↓

5.  Press  Shock at: **Broad complex:** 120j, 150j then 200j  
**Atrial fibrillation:** 200j (up to three times)  
**Other narrow complex:** 75j-120j, 150j then 200j

6.  Remove oxygen, clear and press 

7.  **Confirm SYNC mode**, press **and hold** 

# Adult tachycardia



**UNSTABLE**

**STABLE**  
 Seek expert help

**NO**

**Is the QRS narrow (< 0.12 s)?**

