

TABLE 1: AIRWAY DEVICES

Age (years)	Weight (kg)	Guedel	LMA	iGel	Laryngoscope Blade/CMAC	Air-Q (max size uncuffed ETT)
Neonate	1	3.5	1	1	Mac 0 / Miller 0	0 (3)
	2	5				
	3					
	4					
Infant	5	5.5	1.5	1.5	Mac 1 / Miller 1	1.0 (4.5)
	6					
	7					
	8					
	9					
1 to 6	10	6.5 (for 1-6yr old)	2	2	Mac 1 or 2	1.5 (5)
	11					
	12					
	13					
	14					
	15					
	17					
6 to 10	20	7	2.5	2.5	Mac 2 or 3	2 (5.5)
	25					
	30					
10 to 16	35	8 (12yrs to adult)	3	3	Mac 3	3 (7)
	40					
	45					
	50					
	55					
	60					
			4		Mac 3 or 4	

TABLE 2: ET TUBES: <1y SIZE ON WEIGHT, >1Y SIZE ON AGE

Age (years)	Weight (kg)	ETT Size (Microcuff)	ETT Size (Uncuffed)	Oral ETT Length (cm)
Neonate	1		2 - 2.5	7
	2		2.5 - 3	8
	3	3*	3 - 3.5	9
	4		10	
5	11			
6	11.5			
Infant	7	3.5	3.5-4	12
	8			12.5
	9			13
	10			13.5
	1			4
2	4	4.5	14.5	
3		4.5-5	15	
4	4.5	5	15.5	
5		5-5.5	16	
6	5	5.5	16.5	
7		5.5-6	17	
8	5.5	6	17.5	
9		18		
10	6		18.5	
11			19	
12	6.5		19.5	
13			20	
14	7		20.5	
15				

ALWAYS HAVE ONE LARGER AND ONE SMALLER AVAILABLE

* Microcuff ETT should only be used in term neonates >3kg

TABLE 3: DIFFICULT AIRWAY/ACCESSORY DEVICES: SIZE BASED ON ETT INTERNAL DIAMETER

ETT Size ID	Stylet	Bougie	Suction Catheter	NGT	Fibreoptic Scope		Video Bronchoscope	Cook Airway Exchange Catheter
3	Small	5 Fr	6 Fr	6 Fr	Storz 2.8mm	Olympus 2.2mm	Size 3 (No 162) *SNUG, Remove 15mm Connector*	8 Fr
3.5		10 Fr						
4			Size 3 (No. 162)	11 Fr				
4.5								
5	Medium	15 Fr	10 Fr	10 Fr	Storz 4.0mm	Olympus 4.1mm	Size 4 (No. 164)	14 Fr
5.5								
6			12 Fr	12 Fr				
6.5								
7								
7.5								

TABLE 4: INTUBATING THROUGH 1st GENERATION DISPOSABLE LMA

LMA	Uncuffed ETT	Cuffed ETT	
1	3	LMA under size 3 will not allow passage of pilot balloon	
1.5	3.5		
2	4		
2.5	5		
3	6		5.5