

**Direct Current Cardioversion (DCCV)**  
(See overleaf for algorithm)

1. **Sedation or general analgesia** for conscious patients
2. Attach **3 lead ECG** and defibrillation pads & turn dial to **Manual Defib**

3. Press **Enter Sync** button  and observe sense markers

4. If necessary, press lead  to select lead with most sense markers

5. Press  **Broad complex: 150j, 200j then 360j**  
**Atrial fibrillation: 360j (up to three times)**  
**Other narrow complex: 100j, 200j then 360j**

6. Remove oxygen, clear and press charge button 

7. **Confirm SYNC mode**, Press **and hold** 



# Adult tachycardia

